

## Measuring Titratable Acidity

### What is Total Acidity?

- Total amount of organic acids in a solution (which can include lactic, acetic, tartaric, phosphoric, succinic, citric, etc.)
  - Typically each acid is reported as a concentration in grams per liter (g/L)

### What is Titratable Acidity?

- An approximation of total acidity, measures both associated and dissociated hydrogen ions
  - Measures how much a strong base (ex. sodium hydroxide (NaOH)), it takes to reach a basic pH (typically pH 8.2).
  - This details the total available hydrogen ions and is a more accurate to measure of perceived sourness
  - Is typically reported in either g/L or a percent TA, g/100ml.
  - In beer, this calculation is used to measure lactic acid (the most prominent acid in beer), for other beverages (cider for example) this calculation can be modified to reflect their prominent acid.

### How to Measure:

#### What you will need:

- pH meter
- degassed beer sample
- stir plate with magnetic stir bar
- sodium hydroxide (NaOH) in liquid form (typically sold in 0.1M form)
- Pipettes and glassware, with precision down to 0.1 mL (25 or 50mL buret)
- Gloves

### Procedure:

1. Take the specific gravity of the beer.
2. Take a precise amount of degassed beer (ex. 50ml) in a beaker
3. Using the liquid NaOH and your pH meter, measure a precise amount of the NaOH (usually around 0.1-0.5ml) into your beer sample, stirring each time you add the NaOH
4. Take a pH reading
5. Continue adding the NaOH in 0.1-0.5ml increments until you reach 8.2pH (at this pH you have reached the point where NaOH and lactic acid are equivalent in the solution.)
6. Add up the 0.1-0.5ml increments, or the total volume in ml, it took to get to 8.2pH point. This is your ml 0.1M NaOH number.

### Calculating mL Lactic Acid and g/L of a Specific Acid

- Now we can do math. You will need:
  1. volume of beer
  2. the ml of 0.1M NaOH used to get to 8.2pH
  3. and the specific gravity of your beer.
- There are two different calculations that will provide you ml of lactic acid per 100g of beer (Equation A below) or, another calculation to provide you with g/l as a specific acid (**Equation B, more common**).

#### A. ml Lactic Acid per 100g beer Calculation:

Titratable Acidity (TA) as lactic acid =  
ml 0.1M NaOH x 10/ml of beer x specific gravity

Ex. 50ml of beer at 1.010 specific gravity required 5.6 ml of 0.1M NaOH to reach pH 8.2

$$TA = 5.6\text{ml} \times 10/50 \times 1.010$$

$$TA = 56/50.5 = 1.11$$

Or 1.11ml of 1.0M alkali per 100 g of beer.

#### B. g/l as a Specific Acid Calculation:

(Using the 50ml of beer at 1.010 specific gravity that required 5.6ml of 0.1M NaOH solution)

**1.**

$$\text{Total Acidity (mol/L)} = \text{ml} \times 0.1\text{M} / \text{vol of beer}$$

$$\text{Ex. Total Acidity (mol/L)} = 5.6\text{ml} \times 0.1\text{M} / 50\text{ml} \\ = 0.0112 \text{ mol/L}$$

Now we need to find the g/L of lactic acid:

**2.**

$$\text{g/L Lactic Acid} = \text{TA (mol/L)} \times (90\text{g/mol})^*$$

$$\text{Ex. g/L Lactic Acid} = 0.0112 \text{ mol/L} \times (90\text{g/mol}) \\ = 1.008\text{g/L}$$

\*The 90g/mol is a standard correction for lactic acid

You can also express this number as a % lactic acid by:

**3.**

$$\% \text{ Lactic Acid} = \text{g/L lactic acid} / 1,000\text{g}$$

$$\text{Ex. \% Lactic Acid} = 1.008\text{g/L} / 1,000\text{g} = 0.0018 \\ \times 100 = 0.18\%$$

### References:

ASBC Method of Analysis, Beer Method Number 8

For more information you can visit

<https://www.asbcnet.org/Methods/BeerMethods/Pages/Beer-8-SuppInfo.aspx>